

Automation of the Pillar Biosciences oncoReveal Solid Tumor 22 gene panel (ORST22) on the Biomek NGenius Next Generation Library Prep System

Joseph Goffreda¹, Sean Marris¹, Kirsty McFarland¹, Nicholas Lodato¹, Tyler Buit², ShiPing Zou¹, Zhaohui Wang¹

¹ Pillar Biosciences, Natick, MA; ² Beckman Coulter, Indianapolis, IN



ABSTRACT

Introduction: Fully automated NGS library preparation workflows are highly desirable for clinical settings. In this study, we used the Biomek NGenius instrument (Beckman Coulter) as our platform because it has all necessary components on-deck including cold storage blocks, a magnet position, and a thermal cycler. We developed a fully automated process for our research-use-only (RUO) oncoReveal™ 22-Gene Solid Tumor assay (ORST22) on the Biomek NGenius to prepare 4 to 24 libraries requiring only ~1 hour upfront hands-on time and 6 hours of walkaway time from input DNAs to final libraries. In this study, we evaluated the performance of the automated process and assessed the NGS results against the results from the manual workflow.

Method: Two DNA samples with various replicates were used for the study: NA12878 DNA (genome-in-a-bottle) at ~10ng and the moderately formalin-compromised positive DNA (ModfcDNA, Horizon Discovery) at ~20ng. One to three replicates of no-template control (NTC) were also processed in each run. To evaluate the run-to-run precision, three automation runs were performed and sequenced on different days to prepare 8, 16 and 24 libraries, respectively. To assess the NGS results against those from the manual process, the same reagents and the same input DNAs were used for library preparation by both workflows with 24 libraries for automation and 8 for manual prep. Subsequently, the libraries were sequenced on the same flow cells using Illumina NextSeq 550. The secondary analysis metrics and the variants reported by Pillar's PiVAT bioinformatic software were assessed and compared.

Results: No cross-contamination was observed in any NTC replicates in any of the runs. All runs generated high-quality results with consistent >90% effective on-target rates and ~100% of bases covered at ≥ 0.2x mean coverage. Compared with the manual NGS libraries, the automated libraries showed slightly better uniformity with higher percentages of bases covered at ≥0.4x to 1x of the mean coverage. 20 known variants from each positive ModfcDNA with expected VAF from 1-30% were highly concordant between runs with within runs. No false positive variants were detected in any of the samples.

Conclusion: The automated Pillar ORST22 assay on Biomek NGenius system produced high quality and robust results with minimal hands-on time and human intervention. Together with the easy-to-use setup of the NGenius, and Pillar's single day library preparation workflow, these combined technologies could substantially reduce the lab operation burden.

APPLICATION DESIGN

| | Manual | NGenius |
|--------------------------|--------------------|----------------|
| Sample Count | Up to 24 per batch | 4-24 per batch |
| Input Mass | 10-80ng | |
| Supported DNA Types | Genomic, FFPE | |
| Estimated hands-on time | 3 - 4 hours | 0.5 - 1 hour |
| Estimated hands-off time | 3 hours | 5-6 hours |

Table 1 – Features of ORST22

| App Setting | Description |
|------------------------------------|--|
| Mix beads during Gene-Specific PCR | Mixes AMPure XP beads during gene-specific PCR to prevent them from settling. If not selected, mixing will occur directly before gene-specific product purification. |
| IndexPlate | Allows the operator to enter in a name for the index plate being used in the batch. |
| Indexing PCR Cycles | Allows the operator to set the number of indexing PCR cycles performed within a range of 5-10. |

Table 3 – ORST22 Gene Application settings, configured before each library prep run.

| ORST22 Panel Info | | | | | | | | | | |
|-------------------|--------|------|-------|-------|-------|--------|--------|--------|-------|-------|
| AKT1 | BRAF | DDR2 | ERBB2 | FBXW7 | FGFR2 | KRAS | MET | NRAS | PTEN | STK11 |
| ALK | CTNNB1 | EGFR | ERBB4 | FGFR3 | FGFR3 | MAP2K1 | NOTCH1 | PIK3CA | SMAD4 | TP53 |

Table 2 – The 22 genes covered by the ORST22.

| Gene | Variant ID | Gene | Variant ID |
|------|------------------------------|--------|-----------------------------------|
| EGFR | L858R | PIK3CA | N345I; N345T; N345K |
| | Exon 19 In-Frame Deletions | | E542Q; E542K; E542V |
| | T790M | | E545K; E545Q; E545A; E545G; E545D |
| | G719A; G719C; G719D; G719S | | H1047Y; H1047L; H1047R |
| | Exon 20 In-frame Insertions | | R88Q; R88L |
| KRAS | G12X | FGFR3 | R248C |
| | G13X | | S249C |
| | A59E; A59G; A59T; A59S | | G370C |
| | Q61E; Q61H; Q61K; Q61L; Q61R | | Y373C |
| | K117N | | R130Q; R130L; R130P; R130G; R130* |
| PTEN | A146T; A146P; A146V | PTEN | T319del |
| | E17K | | S249C |
| | BRAF V600E; V600K | | R248C |

Table 4 – Examples of commonly mutated, clinically significant, variants in EGFR, KRAS, PTEN, PIK3CA, FGFR3, AKT1, and BRAF that are targeted in ORST22. Variants in green were used in evaluation during NGenius application development.




Figure 1 – Overview of oncoReveal Solid Tumor 22 Gene Application sections and approximate completion times. For a 24 sample run, the on-deck portion of the workflow can be completed in ~6 hours with several optional stopping points.

RESULTS AND CONCLUSIONS

Figure 2 – Final library yield (nM) across preparation method and sample type conditions.




Figure 3 – Average mapping rate (%) and on-target rate (%) across preparation method and DNA type.




Figure 4 – Average Base coverage depth relative to mean coverage at 0.2x, 0.4x, and 1x across library preparation method and DNA type.




Table 5 – Average variant read frequency % for expected ModfcDNA variants across NGenius and manual sample preparation methods. N=2 (manual), 4(NGenius).

| Gene | HGVS | HGVS | Genomic Position | Expected VAF (%) | Manual | | NGenius | |
|--------|----------------|--------------------|-------------------------|------------------|----------|--------------|----------|--------------|
| | | | | | Hit Rate | VAF (%) | Hit Rate | VAF (%) |
| STK11 | c.816C>T | p.Tyr272= | chr19:1221293-1221293 | 9.0 | 2/2 | 8.01 ± 1.01 | 6/6 | 8.55 ± 1.14 |
| TP53 | c.376-7C>T | N/A | chr17:7578561-7578561 | 26.0 | 2/2 | 29.37 ± 0.71 | 6/6 | 27.42 ± 1.38 |
| KRAS | c.38G>A | p.Gly13Asp | chr12:25398281-25398281 | 15.0 | 2/2 | 14.28 ± 1.51 | 6/6 | 14.78 ± 0.9 |
| KRAS | c.35G>A | p.Gly12Asp | chr12:25398284-25398284 | 6.0 | 2/2 | 6.81 ± 0.52 | 6/6 | 6.94 ± 0.61 |
| ERBB2 | c.2411G>A | p.Gly804Asp | chr17:37881082-37881082 | 5 | 2/2 | 4.55 ± 0.62 | 6/6 | 4.78 ± 0.66 |
| CTNNB1 | c.98C>A | p.Ser33Tyr | chr3:41266101-41266101 | 30 | 2/2 | 32.17 ± 0.51 | 6/6 | 31.01 ± 1.75 |
| CTNNB1 | c.133_135del | p.Ser45del | chr3:41266134-41266136 | 10 | 2/2 | 9.92 ± 1.24 | 6/6 | 9.37 ± 1.38 |
| EGFR | c.1839C>T | p.Ala613= | chr7:55233089-55233089 | 8 | 2/2 | 6.94 ± 0.21 | 6/6 | 7.17 ± 0.47 |
| EGFR | c.2155G>A | p.Gly19Ser | chr7:55241707-55241707 | 24.5 | 2/2 | 24.79 ± 0.41 | 6/6 | 24.43 ± 0.99 |
| EGFR | c.2235_2249del | p.Glu746_Ala750del | chr7:55242465-55242479 | 2 | 2/2 | 1.96 ± 0.06 | 6/6 | 1.86 ± 0.31 |
| EGFR | c.2369C>T | p.Thr790Met | chr7:55249071-55249071 | 1 | 2/2 | 1.28 ± 0.39 | 4/6* | 0.81 ± 0.64 |
| EGFR | c.2573T>G | p.Leu858Arg | chr7:55259515-55259515 | 3 | 2/2 | 3.36 ± 0.93 | 6/6 | 3.85 ± 0.47 |
| MAP2K1 | c.167A>C | p.Gln56Pro | chr15:66727451-66727451 | 30 | 2/2 | 32.29 ± 0.71 | 6/6 | 31.47 ± 2.19 |
| NRAS | c.181C>A | p.Gln61Lys | chr1:1152563011525630 | 12.5 | 2/2 | 11.35 ± 0.08 | 6/6 | 11.51 ± 0.69 |
| BRAF | c.1799T>A | p.Val600Glu | chr7:140453136140453136 | 10.5 | 2/2 | 12.12 ± 1.33 | 6/6 | 13.25 ± 0.95 |
| DDR2 | c.351G>A | p.Met117Ile | chr1:162724579162724579 | 10 | 2/2 | 9.79 ± 1.82 | 6/6 | 9.93 ± 0.92 |
| PIK3CA | c.1173A>G | p.Ile391Met | chr3:178927410178927410 | 9 | 2/2 | 8.56 ± 0.86 | 6/6 | 8.53 ± 0.56 |
| PIK3CA | c.1633G>A | p.Glu54Lys | chr3:178936091178936091 | 9 | 2/2 | 9.11 ± 0.45 | 6/6 | 8.43 ± 0.87 |
| PIK3CA | c.3140A>G | p.His1047Arg | chr3:178952085178952085 | 17.5 | 2/2 | 18.39 ± 0.96 | 6/6 | 17.66 ± 0.75 |
| FGFR3 | c.1953G>A | p.Thr651= | chr4:1807894-1807894 | 100 | 2/2 | 98.83 ± 0.09 | 6/6 | 99.00 ± 0.28 |
| EGFR | c.2361G>A | p.Gln787= | chr7:55249063-55249063 | 15 | 2/2 | 14.13 ± 0.98 | 6/6 | 13.83 ± 1.85 |

*Variant is below VAF cutoff in PiVAT software.

Conclusions:

- Using the NGenius ORST22 application, 24 libraries can be prepared in approximately 7 hours with ~1 hour hands-on time and ~6 hours hands-off time (Figure 1).
- The NGenius automated platform can produce high quality libraries using Pillar Biosciences oncoReveal Solid Tumor 22 Gene panel (Figures 2-4).
- Resulting variant calls for libraries prepared on the NGenius platform were as expected and comparable to manual preparation (Table 5).