

# Thalassemia Panel

The **inheritReveal™ Thalassemia Panel** is specifically designed to cover mutations in the hemoglobin alpha chain (HBA1 and HBA2) and beta chain (HBB), which are responsible for  $\alpha$ -Thalassemia and  $\beta$ -Thalassemia. The panel uses proprietary Stem-Loop Inhibition-Mediated amplification (SLIMamp®) technology, a tiled amplicon-based library prep chemistry for efficient single-tube target enrichment. Additionally, known long deletions (10-30kb) are detected by gap-PCR.

## Panel specifications\*

Enrichment chemistry	Multiplex PCR using tiled amplicons
Number of pools	1 pool
Number of genes/amplicons	3/130
Number of targets	Hotspots in HBA1 and HBA2; Full CDS coverage of HBB (including splice sites), 3' UTRs & 5' UTRs with promoter regions and pathogenic intronic regions; CNAs in HBA1, HBA2, and HBB; 4.8kb total size (excluding gap-PCR amplicons)
Variant types	SNVs, small and medium indels, large deletions (10kb - 30kb), CNAs
Average amplicon size	156bp (excludes gap-PCR amplicons)
Recommended DNA input range	10ng to 75ng
Sample types	Genomic DNA
Mapping rate	99.7% $\pm$ 0.1%
% on-target aligned reads	99.5% $\pm$ 0.1%
Coverage uniformity (% targets with >0.2X mean coverage)	96.5% $\pm$ 2.4%
Recommended reads per sample	~260,000 paired-end reads
Total assay time (from DNA to sequencer)	<7.5 hours

\* Mapping rate, percentage of on-target aligned reads, and coverage uniformity metrics are based on internal testing performed using reference standard materials

### Simple NGS library prep workflow

Maintain control of samples and results with single-tube, tiled amplification that can be performed in-house by any NGS lab


### Sensitive and robust chemistry

Achieve variant detection as low as 1% VAF<sup>†</sup> even with limited DNA input or poor sample quality

### Reduced fully- loaded lab costs

Improve lab efficiency and reduce “no calls”, repeat testing, and difficult interpretation decisions

## Simple, one-day workflow



## Ordering information

Select the panel AND one of the index kit options listed below.

Panel	Part number
inheritReveal™ Thalassemia Panel (24 reactions)	HDA-TH-1001-24

Pillar Index Kit options	Reactions	Part number
Pillar Custom Index Primers Kit A	32 Combinations, 96 reactions	IDX-PI-1001-96
Pillar Custom Index Primers Kit D	96 Combinations, 192 reactions	IDX-PI-1004-192

TO ORDER OR LEARN MORE:  
[pillarbiosci.com](https://pillarbiosci.com)

**For research use only. Not for use in diagnostic procedures.**

© 2022 Pillar Biosciences. Pillar®, SLIMamp®, PiVAT®, oncoReveal® and inheritReveal™ are trademarks of Pillar Biosciences, Inc.  
Current as of 07.2024 | MK-0015-4 v2.0

**PILLAR**  
BIOSCIENCES

Pillar Biosciences, Inc.  
9 Strathmore Rd  
Natick, MA 01760  
(800) 514-9307  
[info@pillarbiosci.com](mailto:info@pillarbiosci.com)